Instructions for use [EN]

LIPOTRIS™

Liposomal intra-articular composite gel, 2 ml

Composition
Each 1 ml contains:
Sodium hyaluronate 22.000 mg
1-palmitoyl-2-oleoyl-phosphatidylcholine (POPC) 0.800 mg
Sodium chloride 9.060 mg
Disodium hydrogen phosphate, anhydrous 0.241 mg
Sodium hydrogen phosphate dihydrate 0.047 mg
Water for injection q.s.

Description
Lipotris™ is a liposomal composite gel for intra-articular application, developed on the medical advice and liposomes. The hyaluronic acid content is obtained by bacterial fermentation and the liposomes are formed from 1-palmitoyl-2-oleoyl-phosphatidylcholine (POPC) of synthetic origin. The two components of the formulation are able to interact with each other to form a composite material - a material with unique tribological properties. The medical device is sterile and non-pyrogenic.

Indications
Pain and decreased articular mobility, associated with traumatic or degenerative changes including osteoarthritis.

Contraindications
The medical device should not be used in patients:
- with known hypersensitivity to any of the components of the product;
- with septic arthritis;
- with dermatosis at the injection site;
- taking anticoagulants.

Precautions
All safety precautions during procedure should be taken to prevent occurrence of septic arthritis.
Medical device intended only for intra-articular injections.
Do not connect it with any disinfectants containing quaternary ammonium salts for skin preparation as hyaluronan can precipitate in their presence.
Do not use the medical device if the sterile barrier is damaged with this medical device or there is no clinical data available on its use in these patients.
Medical device should not be used if a pre-filled syringe is damaged.
Medical device must be used before the expiration date indicated on the packaging.
Do not resterilise.
The pre-filled syringe is intended for single-use only.
Store at temperature between +2°C and +25°C. Keep away from sunlight. Protect from moisture.
Medical device must be kept out of the reach of children.

Potential side effects
Biocompatibility studies have confirmed the high safety of the medical device. Transient local reactions such as pain, redness and swelling may occur after intra-articular injection. These effects can be reduced by cooling the injection site for 5 to 10 minutes after the injection. Simultaneous oral intake of analgesics and anti-inflammatory drugs (NSAIDs) may be beneficial for pain relief.

Serious adverse events reporting
In case of serious adverse events, it should be reported to the manufacturer and to the Member State competent authority. By reporting adverse events more information on safety of the medical device can be gathered.

Interactions
No data.

Administration and dosage
Disinfect the skin at the site of the planned injection. Take the pre-filled syringe from the blister, unscrew the Luer-lock closure and screw a suitable sterile needle onto the syringe. Remove air bubbles before injection. After injection, place the used syringe in a medical waste container and dispose of it.
In case of intra-articular effusion, it is recommended to remove it before the administration of the device.
The therapy consists of a series of three intra-articular injections, administered at weekly intervals. The beneficial effects of the treatment last for at least six months. The treatment cycle may be repeated, if required.

Properties and mode of action
The composite gel Lipotris™ takes advantage of the fact that hyaluronic acid with its precisely selected properties is able to interact closely with liposomes formed from phospholipids and form a uniform liposomal intra-articular composite gel – Lipotris™. When administered into the joint, Lipotris™ activates tribosynovial lubricating layer on the cartilage surface via a dedicated protein. The Lipotris™ mechanism of action involves the transformation of the type of movement between cartilaginous surfaces from sliding to rolling, thereby enabling the creation of surfaces with extremely reduced friction.

Due to its unique mechanism of action, Lipotris™, when administered intra-articularly, substantially alleviates clinical functional deficiency in reducing pain, as well as increased mobility of the joint affected by traumatic or degenerative changes.

Package
Pre-filled syringe with gel placed in a blister, instructions for use, implant card, patient-doctor stickers, sterile ZIG needle.

Needle manufacturer:
Becton Dickinson S.A.
Carretera de Mequinenza, s.n., 22520 Fraga (Huesca), Spain

The needle is sterilized with ethylene oxide.

Symbol	Definition	Symbol	Definition
	Manufacturer	CE	CE marking of conformity
	Keep away from sunlight	LOT	Lot number
	Keep dry		Expiry date
	Temperature limit		Consult instructions for use
	Do not re-sterilise		Caution
	Do not reuse	MD	Medical device
	Do not use if package is damaged and consult instructions for use	UDI	Unique device identifier
STERILE	Sterilized using steam or ethylene oxide		Distributor
	Single sterile barrier system with protective packaging outside		Non-pyrogenic
	Date of manufacture	REF	Reference number
CH REP	Swiss Authorized Representative	UK RP	UK Responsible Person

Not all symbols appear on the labeling.

For more information refer to the Summary of Safety and Clinical Performance available in the Eudamed database.

To be used exclusively by a medical doctor.

Date of issue: 01/2026
No. of the latest revision of the instructions for use: 01/2026

Manufacturer:
Biovico Sp. z o.o.
Hutnicza 15 B
81-061 Gdynia, Poland, EU
+48 58 660 44 88
office@biovico.com
www.biovico.com

Global Distributor:
Implal Sp. z o.o.
Hutnicza 15 B
81-061 Gdynia, Poland, EU
+48 58 727 06 06
office@implal.com
www.implal.com

IMPLAL

INSTRUCTIONS FOR COMPLETING THE IMPLANT CARD [EN]

1. PURPOSE OF THE IMPLANT CARD
The implant card is a document given to the patient after the implantation of a medical device, containing:
- product identification data,
- manufacturer's data,
- data of the implanting person and patient's data.

Its aim is to ensure that the patient and the healthcare system can quickly identify the device and respond effectively to medical incidents.

2. GENERAL PRINCIPLES
• All fields must be completed legibly.
• Product identification information (name, lot number, UDI, expiry date) is affixed by means of a manufacturer's sticker.
• The remaining fields shall be completed by the medical staff after implantation.

3. INSTRUCTIONS – PRODUCT SIDE
To be completed by medical staff after implantation.
• **Device type** – insert "Liposomalny intra-articular composite gel".

Place one of the three stickers supplied with the product containing product identification data in the designated space:
• product name,
• unique device identifier (UDI),
• UDI device identifier (UDI-DI);
• lot number,
• expiry date.

4. INSTRUCTIONS – PATIENT SIDE
To be completed by medical staff after implantation. All fields must be completed legibly, preferably in block capitals.

– enter the patient's full name or patient ID number.

– enter the date of the procedure (DD-MM-YYYY).

– enter the name and address of the facility and the name of the physician performing the implantation.

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